

UPCERA




Hyramic

---

Reliable Materials Expert

# Hyramic

- Available for Single & Multilayer Shades
- Elastic Modulus Approximates to Natural Dentin
- Simple Milling, No Need for Sintering
- Chairside Instant Restoration



Strength  
**≥180 MPa**



## Chemical Composition


Resin polymers	13%-43%
Inorganic filler	55%-85%
Additive	<1%
Colorant	<1%

## Technical Parameters

Water absorption	≤25µg/cm <sup>3</sup>
Solubility in water	≤2µg/cm <sup>3</sup>



## Apply for




Open CAD/CAM System,  
Cerec in lab (Sirona) System

## Indications




Veneers



Inlays / Onlays




Full Contour  
Anterior Crown




Full Contour  
Posterior Crown

## Available for Single & Multilayer Shades

LT 18-14-12



ML 18-14-12



Follow us on:



Upcera\_Official



upcera\_\_official



UPCERA